



National Institute of Health (NIH) Dr. Ricardo Jorge, Lisbon, Portugal



## INSaFLU-TELEVIR: an open web-based bioinformatics suite for metagenomic virus detection and routine genomic surveillance

<https://insaflu.insa.pt>


### Joana Isidro

National Institute of Health Dr. Ricardo Jorge, Portugal  
Department of Infectious Diseases  
Bioinformatics Unit

**[joana.isidro@insa.min-saude.pt](mailto:joana.isidro@insa.min-saude.pt)**

# New era of virus surveillance *Targeting the whole-genome*

Many laboratories do not have bioinformatics capabilities and/or staff needed to analyze raw NGS data, which is an obstacle in the transition to whole-genome sequencing



<https://insaflu.insa.pt>

*Borges et al, 2018*  
*Genome Medicine* 10: 46

# Which are the **main features** of INSaFLU?

<https://insaflu.insa.pt>



**Online platform (free, confidential accounts) but it can be installed locally**



**User-oriented / no advanced expertise in bioinformatics is needed**



**Applicable to NGS data** (single / paired-end data; Illumina / Ion Torrent / ONT)  
**collected from diverse amplicon-based schemas**

# Which are the **main features** of INSaFLU?



<https://insaflu.insa.pt>



Tailored to seasonal human influenza and **SARS-CoV-2** but  
is **applicable to other pathogens**



Dynamically integrates data in a cumulative manner, fitting the  
continuous epidemiological surveillance during epidemics



Outputs are standardized and compatible with other platforms  
for downstream data visualization/analysis



<https://insaflu.insa.pt>

# INSaFLU-TELEVIR: an open web-based bioinformatics suite for metagenomic virus detection and routine genomic surveillance (influenza, SARS-CoV-2 and other pathogens)

## INPUTS


✓ **Sample metadata table**

+

✓ **NGS data \***



*\*Compatible with distinct sequencing technologies:  
Illumina, Ion Torrent and Oxford Nanopore*



<https://insaflu.insa.pt>

## OUTPUTS

### ROUTINE GENOMIC SURVEILLANCE

- ✓ Read quality analysis and improvement
- ✓ Human betacoronaviruses, MPXV and influenza type/subtype classification (directly from reads)
- ✓ Interactive coverage analysis
- ✓ Mutation annotation and consensus sequence generation/curation
- ✓ Gene- and genome-based alignments and phylogenies
- ✓ Intra-host minor variant detection (including flagging of putative mixed infections)
- ✓ SARS-CoV-2 "Pango lineage" classification
- ✓ Direct link to "Nextclade" navigation
- ✓ Genotype-phenotype screening

INSaFLU-TELEVIR development is being co-funded by the European Commission on behalf of OneHealth EJP

TELEVIR project: <https://onehealthejp.eu/jrp-tele-vir/>







Classification

Reads quality analysis and improvement

Mutation detection and consensus generation

Intra-host minor variant detection (and detection of putative mixed infections)

Alignment/Phylogeny

Coverage analysis

SARS-CoV-2 Pango lineage classification and direct link to Nextclade

Integrative Nextstrain phylogenetic and geotemporal analysis.

Screening of Spike aa changes (full repertoire and mutations of interest)

## READ QC

(Trimmomatic / NanoFilt)

- GRAPHICAL QUALITY REPORTS (pre- and post-QC) (FastQC, NanoStat, RabbitQC)
- QUALITY-PROCESSED READS

## REFERENCE-based MAPPING

(BWA / Medaka)

## MUTATION DETECTION AND ANNOTATION

(Freebayes / Medaka, bcftools / SnpEff)

## CONSENSUS GENERATION / CURATION

(automatic masking of low coverage regions; MSA\_masker)

## GRAPHICAL MAPPING PROFILE (IGV)

## COLOR-CODES FOR COVERAGE

## GENOME AND GENE/PROTEIN ALIGNMENTS

(MAFFT, MSASviewer)

## PHYLOGENETIC TREES LINKED TO METADATA VISUALIZATION

(FastTree, PhyloCanvas)

## GENOTYPE-PHENOTYPE SCREENING

## NEXTSTRAIN PHYLOGENETIC AND GEOTEMPORAL ANALYSIS.

(Seasonal Influenza, SARS-CoV-2, Monkeypox, and a “generic” build other viruses)


Influenza A and B type/sub-types  
5 Human BetaCoronaviruses and MPXV  
(SPAdes, ABRicate)


MINOR VARIANTS DETECTION AND ANNOTATION  
(Freebayes, bcftools, SnpEff)


illumina®


PANGO LINEAGE CLASSIFICATION AND NEXTCLADE LINKAGES


→ **Register and Login**  
to a user restricted account


 Home

 References

 Samples

 Projects

 Datasets

 Settings

- ➔ **References:** default reference database; add new references
- ➔ **Samples:** main repository, add/delete samples, inspect / refine QC
- ➔ **Projects:** reference-based mapping and output data navigation
- ➔ **Datasets (new):** Nextstrain phylogenetic and genomic analysis, geotemporal data visualization and exploration of sequence metadata
- ➔ **Settings:** user-defined configuration of parameters for reads quality analysis, mapping, mutation validation and consensus generation/masking

Login, please.

Please, use the browsers Firefox, Safari or Chrome for better view experience.

User name/Email\*

demo

A valid user name or email

Password\*

\*\*\*\*\*

☒ Demo login. Make a demo login. Only can view predefined data.

Login

Register

Reset password

Cancel

[Get another confirmation e-mail \(check your spam first\)](#)

Applicable to:











# INSaFLU projects - reference-based mapping and output data navigation



Instituto Nacional de Saúde

Doutor Ricardo Jorge



INSaFLU

[Documentation](#)

demo\_insaflu - Logout

Home / Projects / Show project results

Project 'project\_SARS\_CoV\_2\_demo'

Phylogenetic trees and metadata visualization

Mutations list

Nucleotide alignments by MSAViewer

Amino acid alignments by MSAViewer

Aln2Pheno report

Intra-host minor variants annotation and uncovering of mixed infections





















Coverage for all samples

Download

Jump to NextClade

Name, Mixed, Type and D 


Search

Sample Name	Type And Subtype	Putative Mixed Infection	Technology	Dataset	Coverage	Consensus File	Alerts	Results
 demo_SARSCoV2_001	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2	<div><div></div></div>	<a href="#">Consens...01.fasta</a>	1	<div><div></div><div></div><div> More info</div></div>
 demo_SARSCoV2_002	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2	<div><div></div></div>	<a href="#">Consens...02.fasta</a>	1	<div><div></div><div></div><div> More info</div></div>
 demo_SARSCoV2_003	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2	<div><div></div></div>	<a href="#">Consens...03.fasta</a>	1	<div><div></div><div></div><div> More info</div></div>
 demo_SARSCoV2_004	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2	<div><div></div></div>	<a href="#">Consens...04.fasta</a>	1	<div><div></div><div></div><div> More info</div></div>
 demo_SARSCoV2_005	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2	<div><div></div></div>	<a href="#">Consens...05.fasta</a>	1	<div><div></div><div></div><div> More info</div></div>


Multiple outputs in nomenclature-stable and standardized formats that can be downloaded or visually explored *in situ*.



# Examples of INSaFLU projects outputs



Instituto Nacional de Saúde  
Doutor Ricardo Jorge



INSaFLU

Documentationdemo\_insaflu - Logout

Home / Projects / Show project results

Project 'project\_SARS\_CoV\_2\_demo'

Phylogenetic trees and metadata visualization

Mutations list

Nucleotide alignments by MSAViewer

Amino acid alignments by MSAViewer


























Aln2Pheno report

Intra-host minor variants annotation and uncovering of mixed infections

Coverage for all samples

DownloadJump to NextClade

Name, Mixed, Type and DSearch

Sample Name	Type And Subtype	Putative Mixed Infection	Technology	Dataset	Coverage	Consensus File	Alerts	Results
 demo_SARSCoV2_001	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2		<a href="#">Consens...01.fasta</a>	1	   More info
 demo_SARSCoV2_002	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2		<a href="#">Consens...02.fasta</a>	1	   More info
 demo_SARSCoV2_003	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2		<a href="#">Consens...03.fasta</a>	1	   More info
 demo_SARSCoV2_004	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2		<a href="#">Consens...04.fasta</a>	1	   More info
 demo_SARSCoV2_005	BetaCoV-SARS_CoV_2	No	Illumina/IonTorrent	demo_sarscov2		<a href="#">Consens...05.fasta</a>	1	   More info

Expand-and-collapse panels

Multiple outputs in nomenclature-stable and standardized formats that can be downloaded or visually explored *in situ*.

# Examples of “Projects” outputs



<https://insaflu.insa.pt>

Mutations in consensus




SAMPLE ID	CHROM	POS	TYPE	REF	ALT	EVIDENCE	FTYPE	STRAN	NT_PC	AA_PC	EFFECT	NT CHANGE	AA CHANGE	LOCUS	GENE	PRODUCT
SARS_CoV_2_demo_01	MN908947	241	snp	C	T	T:5237 C:5	5'UTR	+								
SARS_CoV_2_demo_01	MN908947	3037	snp	C	T	T:5634 C:3	CDS	+	2772/132 924/4400		synonymous_variant	2772C>T	Phe924Phe		orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_01	MN908947	14408	snp	C	T	T:2898 C:1	CDS	+	941/8088 314/2695		missense_variant	941C>T	Pro314Leu		orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_01	MN908947	23403	snp	A	G	G:10357 A:2	CDS	+	1841/382 614/1273		missense_variant	1841A>G	Asp614Gly		S	surface glycoprotein
SARS_CoV_2_demo_01	MN908947	608	snp	GGG	AAC	AAC:3205 GGC	CDS	+	608/1260 203/419		missense_variant	608_610delIGGGinsAAC	ArgGly203LysArg		N	nucleocapsid phosphoprotein
SARS_CoV_2_demo_01	MN908947	3461	snp	C	T	T:3461 C:3	5'UTR	+								
SARS_CoV_2_demo_01	MN908947	4320	snp	C	T	T:4320 C:4	CDS	+	2772/132 924/4400		synonymous_variant	2772C>T	Phe924Phe		orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_01	MN908947	4503	snp	C	T		CDS	+	4278/212 1426/709		synonymous_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_01	MN908947	10223	snp	C	T		CDS	+	9998/212 3333/709		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_01	MN908947	25284	snp	C	A		CDS	+	3762/382 1254/127		stop_gained				S	surface glycoprotein
SARS_CoV_2_demo_01	MN908947	28999	snp	A	T		CDS	+	766/1260 256/419		stop_gained				N	nucleocapsid phosphoprotein
SARS_CoV_2_demo_02	MN908947	2863	snp	A	G		CDS	+	2638/212 880/7095		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_02	MN908947	8168	snp	C	T		CDS	+	7943/212 2648/709		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_02	MN908947	24328	snp	G	T		CDS	+	2806/382 936/1273		missense_variant				S	surface glycoprotein
SARS_CoV_2_demo_02	MN908947	25284	snp	C	A		CDS	+	3762/382 1254/127		stop_gained				S	surface glycoprotein
SARS_CoV_2_demo_02	MN908947	28999	snp	A	T		CDS	+	766/1260 256/419		stop_gained				N	nucleocapsid phosphoprotein
SARS_CoV_2_demo_03	MN908947	3907	snp	G	T		CDS	+	3682/212 1228/709		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_03	MN908947	3926	snp	C	T		CDS	+	3701/212 1234/709		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_03	MN908947	25284	snp	C	A		CDS	+	3762/382 1254/127		stop_gained				S	surface glycoprotein
SARS_CoV_2_demo_03	MN908947	25310	snp	C	T		CDS	+	3788/382 1263/127		missense_variant				S	surface glycoprotein
SARS_CoV_2_demo_04	MN908947	1857	snp	A	G		CDS	+	1632/212 544/7095		synonymous_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_04	MN908947	24321	snp	A	G		CDS	+	2799/382 933/1273		synonymous_variant				S	surface glycoprotein
SARS_CoV_2_demo_04	MN908947	25284	snp	C	A		CDS	+	3762/382 1254/127		stop_gained				S	surface glycoprotein
SARS_CoV_2_demo_04	MN908947	25476	snp	C	T		CDS	+	124/828 42/275		missense_variant				ORF3a	ORF3a protein
SARS_CoV_2_demo_04	MN908947	28999	snp	A	T		CDS	+	766/1260 256/419		stop_gained				N	nucleocapsid phosphoprotein
SARS_CoV_2_demo_05	MN908947	64	snp	T	C											
SARS_CoV_2_demo_05	MN908947	789	snp	C	T		CDS	+	564/2129 188/7095		synonymous_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_05	MN908947	6571	snp	G	A		CDS	+	6346/212 2116/709		missense_variant				orf1ab	orf1ab polyprotein
SARS_CoV_2_demo_05	MN908947	25006	snp	C	T		CDS	+	3484/382 1162/127		missense_variant				S	surface glycoprotein
SARS_CoV_2_demo_05	MN908947	25284	snp	C	A		CDS	+	3762/382 1254/127		stop_gained				S	surface glycoprotein
SARS_CoV_2_demo_05	MN908947	28999	snp	A	T		CDS	+	766/1260 256/419		stop_gained				N	nucleocapsid phosphoprotein
SARS_CoV_2_demo_06	MN908947	2073	snp	C	T		CDS	+	1848/212 616/7095		synonymous_variant				orf1ab	orf1ab polyprotein


minor variants

## Examples of “Projects” outputs

# Genome alignment



## Protein alignments






e.g. ORF1ab

e.g. Spike



# Examples of “Projects” outputs

## Phylogenetic trees and metadata visualization



# Examples of “Projects” outputs

*align2pheno* tool (pilot module for SARS-CoV-2) screens and reports mutations of interest in Spike.

Aln2Pheno report					
Sequence	lineage	Flagged mutations	All mutations	Nflagged	Nmutations
	Clear	Clear			
demo_sarscov2_001	B.1.525	S:H69del;S:V70del;S:Y144del;S:E484K	S:Q52R;S:A67V;S:H69del;S:V70del;S:Y144del;S:E484K;S:A570V;S:D614G;S:Q677H;S:F888L	4	10
demo_sarscov2_002	AY.102	S:T95I;S:R158G;S:L452R	S:T19R;S:T95I;S:E156del;S:F157del;S:R158G;S:L452R;S:T478K;S:D614G;S:P681R;S:D950N	3	10
demo_sarscov2_003	AY.1	S:T95I;S:R158G;S:K417N;S:L452R	S:T19R;S:T95I;S:E156del;S:F157del;S:R158G;S:W258L;S:K417N;S:L452R;S:T478K;S:D614G;S:P681R;S:D950N	4	12
demo_sarscov2_004	P.1	S:L18F;S:K417T;S:E484K;S:N501Y;S:H655Y	S:L18F;S:T20N;S:P265;S:D138Y;S:R190S;S:K417T;S:E484K;S:N501Y;S:D614G;S:H655Y;S:T1027I;S:V1176F	5	12
demo_sarscov2_005	B.1.351.3	S:L18F;S:K417N;S:E484K;S:N501Y	S:L18F;S:A67V;S:D80A;S:D215G;S:L241del;S:L242del;S:A243del;S:K417N;S:E484K;S:N501Y;S:D614G;S:A701V	4	12
demo_sarscov2_012	B.1.1.7	S:H69del;S:V70del;S:Y144del;S:N501Y	S:H69del;S:V70del;S:Y144del;S:N501Y;S:A570D;S:D614G;S:P681H;S:T716I;S:S982A;S:D1118H	4	10
demo_sarscov2_013	B.1.525	S:H69del;S:V70del;S:Y144del;S:E484K	S:Q52R;S:A67V;S:H69del;S:V70del;S:Y144del;S:E484K;S:D614G;S:Q677H;S:F888L	4	9
demo_sarscov2_014	B.1.1.7	S:H69del;S:V70del;S:Y144del;S:N501Y	S:P25L;S:H69del;S:V70del;S:Y144del;S:N501Y;S:A570D;S:D614G;S:P681H;S:T716I;S:S982A;S:D1118H	4	11
demo_sarscov2_016	AY.124	S:T95I;S:R158G;S:L452R	S:T19R;S:T95I;S:E156del;S:F157del;S:R158G;S:L452R;S:T478K;S:D614G;S:P681R;S:D950N	3	10
demo_sarscov2_017	B.1.525	S:H69del;S:V70del;S:Y144del;S:E484K	S:Q52R;S:A67V;S:H69del;S:V70del;S:Y144del;S:E484K;S:A570V;S:D614G;S:Q677H;S:F888L	4	10


Screening of SARS-CoV-2 Spike amino acid alignments against two “genotype-phenotype” databases: the COG-UK Antigenic mutations (<https://sars2.cvr.gla.ac.uk/cog-uk/>) and the Pokay Database (<https://github.com/nodrogluap/pokay/tree/master/data>).





<https://cov-lineages.org/>

## Automate (and up-to-date) Pango lineage classification




## Direct link to Nextclade navigation




Nextclade<sub>v1.8.1</sub>  
Clade assignment, mutation calling, and sequence quality checks

<https://clades.nextstrain.org/>



# DATASETS – New module



**References:** default reference database; add new references


**Samples:** main repository, add/delete samples, inspect / refine QC

**Projects:** reference-based mapping and output data navigation

**Datasets (new):** Nextstrain phylogenetic and genomic analysis, geotemporal data visualization and exploration of sequence metadata

**Settings:** user-defined configuration of parameters for reads quality analysis, mapping, mutation validation and consensus generation/masking

→ **Register and Login**  
to a user restricted account



Login, please.

Please, use the browsers Firefox, Safari or Chrome for better view experience.

User name/Email\*

demo

A valid user name or email

Password\*

\*\*\*\*\*

☒ Demo login. Make a demo login. Only can view predefined data.

Login Register Reset password Cancel

[Get another confirmation e-mail \(check your spam first\)](#)

Applicable to:

illumina®

iontorrent

Oxford  
**NANOPORE**  
Technologies



[Home](#) / [Datasets](#)
[+ New Dataset](#)



Dataset Name	NextStrain Build	Last Change Date	Creation Date	Alerts	#Sequences (P/C/R)	Results
<a href="#">Nextstrain_MPXV_demo</a>	Monkeypox	10-11-2022 23:02	10-11-2022 22:47	0	<a href="#">(7/0/1)</a> <a href="#">Add Sequences</a>	<a href="#">i</a> <a href="#">p</a> <a href="#">a</a>
<a href="#">Nextstrain_SARS_CoV_2_demo</a>	SARS-CoV-2	10-11-2022 20:00	09-11-2022 22:15	0	<a href="#">(66/0/0)</a> <a href="#">Add Sequences</a>	<a href="#">i</a> <a href="#">p</a> <a href="#">a</a>
<a href="#">Nextstrain_influenza_A_H3N2_demo</a>	Influenza (H3N2 HA 12years)	10-11-2022 20:04	09-11-2022 22:14	0	<a href="#">(13/0/0)</a> <a href="#">Add Sequences</a>	<a href="#">i</a> <a href="#">p</a> <a href="#">a</a>
<a href="#">Nextstrain_generic_RSV_demo</a>	Generic	10-11-2022 20:05	09-11-2022 22:13	0	<a href="#">(4/0/1)</a> <a href="#">Add Sequences</a>	<a href="#">i</a> <a href="#">p</a> <a href="#">a</a>

[Refresh](#)

Total datasets: 4

Available builds: Seasonal Influenza, SARS-CoV-2, Monkeypox and Generic

# Examples of “Datasets” outputs



<https://insaflu.insa.pt>

Nextstrain (<https://nextstrain.org/>) phylogenetic and geotemporal analysis.

Powered by




Dataset 'SARS\_CoV\_2\_demo'

Phylogenetic trees and metadata visualization

Nucleotide alignments by MSAViewer



Interactive exploration of Nextstrain alignments and phylogenetic trees using PhyloCanvas



Download Metadata for NextStrain

Download - Dataset\_list.csv

Download - Dataset\_list.tsv

Download - AllConsensus.fasta

Download - Nextstrain\_metadata.tsv


Download - auspice\_json.zip

Download - AllFiles.zip



Integrative and interactive Nextstrain phylogenetic and geotemporal analysis.

DRAG&DROP onto <https://auspice.us/>






# Examples of “Datasets” outputs



<https://insaflu.insa.pt>

Nextstrain (<https://nextstrain.org/>) phylogenetic and geotemporal analysis.

Powered by



*Nextstrain builds for seasonal Influenza, SARS-CoV-2 and Monkeypox as well as a “generic” build for other viruses.*



## INSaFLU-TELEVIR has been continuously updated

- Firstly tailored for **influenza surveillance**
- **SARS-CoV-2** – adapted at beginning of pandemic
- **2022 Monkeypox virus outbreak**

Virus-specific routine  
genomic surveillance



## New challenges - viral metagenomics

- **Virus detection in complex biological samples**
- **Crucial for early virus detection** (especially, novel threats)
- Bioinformatics analysis is a major obstacle/need in the field

Virus detection in  
complex samples



SOURCE: <https://www.insideprecisionmedicine.com/magazine-editions/volume-6-issue-number-4-july-august-2019/mining-the-microbiome-microbiome-startups-are-proliferating-exploring-new-markets-and-colonizing-diverse-therapeutic-niches/>

# TELEVIR Module - Project TELEVIR (OneHealth EJP)

- Implementation of a new module for virus detection from reads

The screenshot displays the INSaFLU web application interface. At the top left is the INSaFLU logo. Below it is a sidebar with navigation icons. The main content area shows a breadcrumb trail 'Home / Project Index' and a section titled 'Modules'. There are three module cards:

- Virus Detection (TELEVIR)**: Features an orange magnifying glass icon. Description: 'From Reads to viral detection'. A blue 'Projects' button is at the bottom.
- Reference-based genome assembly (INSaFLU)**: Features a blue circular icon with horizontal lines. Description: 'From reads to consensus sequence, mutations list, alignments, Pango lineage classification, etc.'. A blue 'Projects' button is at the bottom.
- Nextstrain**: Features a green circular icon with a colorful globe. Description: 'From consensus sequences to phylogeographic and temporal analysis and metadata navigation'. A blue 'Datasets' button is at the bottom.

At the bottom, two brackets group the modules:

- An orange bracket under the 'Virus Detection' module is labeled 'New virus detection module'.
- A grey bracket under the 'Reference-based genome assembly' and 'Nextstrain' modules is labeled 'Genome-based surveillance-oriented modules'.



<https://insaflu.insa.pt>

# INSaFLU-TELEVIR: an open web-based bioinformatics suite for metagenomic virus detection and routine genomic surveillance (influenza, SARS-CoV-2 and other pathogens)

New virus detection module

SOON !

## INPUTS


✓ **Sample metadata table**

+

✓ **NGS data \***



*\*Compatible with distinct sequencing technologies:  
Illumina, Ion Torrent and Oxford Nanopore*



## OUTPUTS

### VIRUS DETECTION (TELEVIR module)

- ✓ Read quality analysis and improvement
- ✓ **State of the art software**
- ✓ Modular Pipeline
- ✓ **Multiple Classification methods**
- ✓ **Multiple Viral databases**
- ✓ **Summary statistics and intuitive and interactive end-user reports**

INSaFLU-TELEVIR development is being co-funded by the European Commission on behalf of OneHealth EJP


TELE-Vir project: <https://onehealthejp.eu/jrp-tele-vir/>



# TELEVIR Module - Project TELEVIR (OneHealth EJP)

## Virus detection from reads

Remove reads from host



Keep the viral reads in the sample

Assembly; taxonomic classification of reads and contigs

Reads/contigs are remapped against reference genomes of the identified taxa

Summary statistics of the identified hits and intuitive end-user reports

# TELEVIR Module - Project TELEVIR (OneHealth EJP)

## *Virus detection from reads*

Run

Pre-Processing

Assembly

Reads and Contigs Classification

Remapping

Pathogen Identification

Multiple useful statistics in report

[Download as CSV](#)

Description	Taxid	accID	Cov (%)	Depth	DepthC	Mapped reads	start prop (%)	mapped_prop (%)	Gaps	Windows Covered	class. success	mapping success	Warning
Influenza A virus (A/Korea/426/1968(H2N2)) segment 8, complete sequence	488241	NC_007380.1	98.57	10.17	11.42	16	0.15035	0.48222	4	3/3	reads	reads	
Influenza A virus (A/New York/392/2004(H3N2)) segment 1, complete sequence	335341	NC_007373.1	98.42	46.89	47.73	87	0.81752	2.62206	2	3/3	reads and contigs	reads	
Influenza A virus (A/Korea/426/1968(H2N2)) segment 1, complete sequence	488241	NC_007378.1	98.33	96.01	97.81	188	1.76659	5.66606	3	3/3	reads	reads	
Influenza A virus (A/New York/392/2004(H3N2)) segment 2, complete sequence	335341	NC_007372.1	97.57	418.13	439.14	963	9.04905	29.02351	2	3/3	reads and contigs	reads	

Horizontal coverage


Depth of coverage

N and % of reads mapped



# TELEVIR Module - Project TELEVIR (OneHealth EJP)

## *Virus detection from reads*




Each hit is an “expand-and-collapse” panel that can be further visualized and explored

- Visualization of coverage throughout reference genome
- Mapped reads and contigs can be downloaded for further confirmations

# TELEVIR Module - Project TELEVIR (OneHealth EJP)

## *Virus detection from reads*



For further information from viral genome (if there is sufficient coverage), the user can **analyze the sample with the “surveillance-oriented” modules in INSaFLU-TELEVIR** (reference-based mapping, phylogenetic trees, ...)

# Acknowledgements - INSaFLU-TELEVIR platform development

Instituto Nacional de Saúde  
Doutor Ricardo Jorge



## Bioinformatics Unit

**Vítor Borges**

**Daniel Sobral**

**João Santos**

Miguel Pinto

Rita Ferreira

Luis Coelho

Verónica Mixão

Ricardo Pais – former Pos-Doc

Alexandra Nunes

**João Paulo Gomes**

## External partners

Miguel Pinheiro (iBiMED / Universidade de Aveiro)

(co-developer)

## OneHealth EJP TELE-Vir project partners

Infraestrutura Nacional de Computação Distribuída (INCD),  
for providing computational resources for testing

## National Reference Laboratory (NRL) for Influenza

### Raquel Guiomar and collaborators

INSaFLU development is being co-funded by the European Commission on behalf of OneHealth EJP

TELE-Vir project: <https://onehealthejp.eu/jrp-tele-vir/>

